



Transatlantic

Produkte für eine heile Welt

Multi- Umbettungstuch

The absorbent pad, carrying sheet and transfer aid with the practical handles.



Transatlantic

Produkte für eine heile Welt

Transatlantic
Handels-Gesellschaft
Stolpe & Co. mbH

Siemensstraße 7
61267 Neu-Anspach



Telefon: +49 6081 943050
Fax: +49 6081 943080

Mail: info@transat.de
www.transatlantic.de

Transatlantic Multi-Umbettungstuch

The revolutionary Multi sheet is not only used for lifting, carrying and pulling patients weighing up to 200 kg, but is also ideally suited as a waterproof pad (with a capacity of approx. 4 liters) and for back-friendly transfers.

Example of use:



1 John had an accident and is bleeding heavily. Immediately, the paramedics put John on the Multi sheet and then on the stretcher.

2 In the emergency room, John is transferred from the stretcher to the hospital bed with the help of a rollboard and the Multi sheet.



3 The doctors take care of John and his open wounds. During the entire treatment, the Multi sheet remains underneath him.

4 John has to get surgery. Only now, the Multi sheet is removed, as it must be sterile and clean in the OR.

✓ Benefits:

- ✓ Absorbent pad and transfer aid with handles in one product
- ✓ Carrying capacity of up to 200 kg
- ✓ Holds up to 4 litres of liquid
- ✓ Safe and stable thanks to 4 handles on both sides
- ✓ Can be used alone or together with our Transaroll boards
- ✓ Fast and clean use
- ✓ Single-use product

Multi-Umbettungstuch

80 x 210 cm, 20 pcs./box
Article No. 4700181K

