

Transaroll®

Patient Transfer

Practical Aids for Easy and Ergonomic Patient Transfer

- ✓ Roll boards & accessories
- ✓ Covers
- ✓ Storage systems
- ✓ Sheets
- ✓ Transfer boards & accessories

 **Transatlantic**
Products for a healthier world



Rollboards & accessories

4

Aids for lateral transfer



Covers

9

Protective covers for hygienic and easy use



Storage systems

10

Hygienic storage and easy transport of the rollboards



Sheets

11

Absorbent transfer aids



Transfer boards & accessories

12

Transfer aids for the wheelchair

Rollboards – aids for lateral transfer

	Silver	Hygiene	Blue	Professional	Smart	Klipp-Klapp
Sizes	5	3	4	5	2	1
Ends	Integrated carrying handles	Integrated carrying handles	Integrated carrying handles	Plastic end strip with loops	End strip with loop	none
Load	350 kg	300 kg	300 kg	200 kg	200 kg	180 kg
Radiolucent	✓	✓	✓	✓	✓	✓
Durability	★★★★★	★★★★★	★★★	★★★	★★	★★
Gliding ability	★★★★★	★★★	★★★★★	★★★★★	★★★★	★
Possible applications	★★★★★	★★★	★★★★★	★★★★★	★★★★	★
Hygiene	★★	★★★★★	★★★	★★	★★	★★★★★
	#1	#2	#3	#4	#5	#6

For more information on the available variants, please refer to the individual product information or our website at: www.transatlantic.de



What does this mean? Standard, Mini, Ambu, Long, Special and Slim are the size units of our rollboards.

**Please note that the sizes of all rollboards may vary by +/- 3%.*

Standard: ca. 180 x 50 cm.

Our most common size, since the boards can be used in all areas.

Mini: between 95-82 x 50 cm.

The boards are primarily used in gynecology and proctology, as they allow transfer even when the legs are elevated.

Ambu: between 86-75 x 45-40 cm.

Because of their narrow size, they are particularly suitable for short distances and transfers such as from narrow stretchers in emergency services. In addition, they are very space-saving when stored and can also be stored under a stretcher.

Long: approx. 195 x 50 cm.

Particularly suitable for large patients and for beds without disturbing the headboard and footboard.

Special: approx. 150 x 60 cm.

Particularly suitable for obese patients and transfers where a longer distance must be crossed.

Slim: approx. 120 x 40 cm.

Particularly suitable for stretcher transfers during emergency room admission and in CT.

For a smooth patient transfer, the heaviest points of the body—shoulders and pelvis—must be on the rollboard. That's why our Transaroll boards are most common in the Standard size. They can be used in all areas and allow for a comfortable transfer for both patient and user, since both the head and the heels rest on the board. But our portfolio also provides special solutions for the operating room with leg elevation (Mini), for the outpatient clinic and emergency room (Ambu, Slim), and for very tall or obese patients (Long, Special).

Basically, however, it is the width of the rollboard that determines its suitability for use. The wider the rollboard, the greater the distances that can be covered. Whether a rollboard is foldable or not has no influence on its transfer capabilities. They can only be stored in a more space-saving way, such as in our wall mount organizer.



Transaroll Silverboard* - the rollboard with the silver cover

The Transaroll Silverboard is the new standard for all lateral positioning and transfer applications. Due to its extremely stable core, the Silverboard allows for ergonomic and easy patient transfer.

Standard (foldable)
180 x 50 cm
Item no. 4500030

Special (not foldable)
150 x 60 cm
Item no. 4500034

Mini (not foldable)
95 x 50 cm
Item no. 4500029

Special (foldable)
150 x 60 cm
Item no. 4500033

Mini (foldable)
95 x 50 cm
Item no. 4500028

Key facts

 **Core**
Strong polyethylene core

 **Cover**
Hygienic, antibacterial gliding fabric

 **Ends**
Integrated carrying handles, extremely stable

 **Load**
350 kg max.

Overview of the benefits:

- ✓ Ideal for obese patients
- ✓ The only rollboard which makes it possible to lift a person
- ✓ No sharp edges, all corners are rounded off
- ✓ No magnetic parts, plastic parts or screws
- ✓ Radiolucent and suitable for MRI
- ✓ Transfer even with soft positioning mattresses
- ✓ No end strip (maximum hygienic safety)
- ✓ Can be disinfected using all common disinfection agents

Transaroll Hygieneboard - maximum hygienic safety

The Hygieneboard ensures the most hygienic use of all. The rollboard consists of a stable and easy-to-clean core and our green gliding cover, which is applied to the rollboard each time it is used.



Standard (foldable)
180 x 50 cm
Item no. 4500026

Mini (foldable)
95 x 50 cm
Item no. 4500036

Mini (non-foldable)
95 x 50 cm
Item no. 4500035

Key facts



Core

Strong polyethylene core



Cover

Disposable gliding cover



Ends

Integrated carrying handles, extremely stable



Load

300 kg max.



Transaroll Blueboard - the rollboard for care

Adapted to the average size and weight of patients in Europe, it is suitable for low-impact transfer in home care, as well as for use in clinics, hospitals and nursing homes that rely on quality and experience. The cover of the rollboard is machine washable.

Standard (foldable)
180 x 50 cm
Item no. 4501100

Long (foldable)
195 x 50 cm
Item no. 4500406

Mini (not foldable)
95 x 50 cm
Item no. 4500400

Mini (foldable)
95 x 50 cm
Item no. 4500403

Key facts



Core

Strong polyethylene core



Cover

Ripstop nylon




Ends

Integrated carrying handles, extremely stable



Load

300 kg max.



Transaroll Professional board* – the classic product for more than 30 years

Our Transaroll Professional boards have been on the German medical market for over 30 years and are as reliable today as they were then. They consist of a special lightweight foam, which is also very durable and resilient despite its low weight.

Standard (not foldable)
180 x 50 cm
Item no. 4401150

Standard (foldable)
180 x 50 cm
Item no. 4401100

Mini (not foldable)
82 x 50 cm
Item no. 4400400

Ambu (not foldable)
80 x 40 cm
Item no. 4400114

Slim (not foldable)
120 x 40 cm
Item no. 44004200

Key facts



Core

Special lightweight foam



Cover

Hygienic, antibacterial gliding sheet



Ends

Plastic endstrip with loops



Load

200 kg max.

Transaroll Smartboard* - the low-cost alternative

This rollboard is perfect for those who are looking for an inexpensive rollboard, but still want to enjoy high gliding capability.

Standard (foldable)
180 x 50 cm
Item no. 4400102

Ambu (foldable)
86 x 40 cm
Item no. 4400101



Key facts



Core

Plate welded in polyester fabric



Cover

Hybrid nylon cover



Ends

End strip with loop



Load

200 kg max.

Transaroll Klipp-Klapp Board - the roll board for every ambulance



Transaroll Klipp-Klapp Board

The Transaroll Klipp-Klapp Board was developed with emergency services in mind. Due to the double foldable body of the board, it allows for quick and easy space-saving storage in the ambulance.

Ambu (double foldable)
75 x 45 cm
Item no. 4400103

Key facts

 **Core**
HDPE Board

 **Cover**
Grey gliding fabric

 **Ends**
No grip bars

 **Load**
max. 180 kg



For emergency services:

Emergency services carry out the initial care in emergency cases. Therefore, the procedures must flow in a rapid, hygienic and uncomplicated manner. Our products are ideally suited for this purpose since they are easy to handle and meet the highest hygiene and quality standards.

Covers - for fast, hygienic use

Our disposable covers and gliding covers are crafted for seamless patient transfers and enhanced hygiene. Designed to protect rollboards from contamination, these covers ensure quick, hygienic, and efficient handling in hospitals and emergency settings. The versatile gliding cover stands out for its ability to assist in patient movement independently, offering an all-in-one solution for diverse healthcare environments.

Overview of the benefits:

- ✓ Gliding cover enables patient transfers without a rollboard
- ✓ Disposable covers protect rollboards from contamination
- ✓ Perforated rolls and folds allow customizable sizing
- ✓ Anti-static options ensure safe handling
- ✓ Single-use design maintains hygiene standards
- ✓ Compact and versatile for easy storage and use

Disposable Gliding Cover

Due to the exceptional gliding properties, the cover is suitable for moving the patient in bed even without a rollboard.

Standard, 60 pcs.

(Roll, perforated)

180 x 53 cm

Item no. 4400204

Disposable, by the meter

(cardboard)

100 m x 60 cm

Item no. 4400203

Disposable, by the meter

Roll, can be cut to size)

100 m x 60 cm

Item no. 4400200



Disposable Protective Cover

Protects the rollboard from possible contamination in individual cases.

Disposable protective cover, 150 pc, anti-static

(Roll, perforated)

54 x 100 cm

Item no. 4500126

Disposable protective cover, 50 pc, standard

(Cardboard, folded)

60 cm x 195 cm

Item no. 4500120

Disposable protective cover, 100 pc, anti-static

(Roll, perforated)

56 cm x 195 cm

Item no. 4500125N

Hygienic storage solutions for safe and clean storage of rollboards and covers

Our wall mounts offer innovative solutions for storing rollboards and protective covers in clinics and care facilities. Whether permanently mounted on the wall or flexible with a mobile suction hook - our products enable hygienic and ready-to-hand storage.

They have been specially developed to meet the requirements of the medical environment and make daily work easier.

Mobile suction hook for hanging up rollboards

Item no.: 4400086



Dispenser for disposable protective cover and suspension for rollboards

Item no.: 4400088



Overview of the benefits:

- ✓ Safe and clean storage of rollboards and protective covers.
- ✓ Space-saving storage.
- ✓ Simple installation of the wall bracket with or without screws.
- ✓ Flexible use thanks to mobile suction hooks for smooth surfaces.
- ✓ Quick access to rollboards and covers, when they are needed

Sheets - absorbent transfer aids

The revolutionary sheets are not only used for lifting, carrying and extricating patients, but are also ideal as waterproof padding and for transferring a patient in a way that protects the back.

Overview of the benefits:

- ✓ Padding and transfer cloth with handles in one product
- ✓ Safe and stable gripping due to multiple handles
- ✓ Can be used alone or together with our Transaroll boards
- ✓ Quick and clean use
- ✓ Disposable product

Transatlantic Multisheet

80 x 210 cm

Item no. 4700181K

Key facts

Material

Outer surface:
Polyethylene cover,
blue

Absorbent core:
Superabsorbent, fluff-
tissue mixture

Fleece: Polypropylene
(water absorbent or
water permeable)

Grip holes


8 integrated grip holes

Fluid uptake

4 l fluid uptake

Load

200 kg max.



XXL Carry Sheet

107 x 210 cm

Item no. 4700081K

Key facts

Material

Polyester
antimicrobial fabric

Fluid uptake

4 l fluid uptake

Grip holes

16 integrated grips

Load

350 kg max.

Transfer boards & accessories - transfer aids for the wheelchair

Our broad portfolio also offers useful aids for seated transfers. Our glide boards are ideally suited for all care areas in clinical and nursing care applications and facilitate ergonomic working.



Silver glider

Classic glider for all applications



Beasy glide board

Nursing glide board with turntable disc



meinHOLTZ transfer board

Design product with individual motifs

What is a transfer board?

Transfer boards make it easier to transfer mobility-impaired people. They form a bridge between the wheelchair and another surface, such as the bed, toilet or car seat.

Who can use a transfer board?

In principle, anyone who finds it difficult to transfer from a wheelchair can use a glide board. It supports an active transfer with the greatest possible independence. A glide board can also be a great support for relatives or people who help people in need of care, and can relieve strain on the back.



Silver Glider - the classic glider in all areas

The M and L Silver Gliders have a shape developed by therapists, which has proven its worth for many years. Due to a very smooth surface, the gliders minimize frictional resistance while anti-slip pads prevent the board from slipping away on the base.

M Silver Glider

67 x 22 cm

Item no. 4500414

L Silver Glider

85 x 25 cm

Item no. 4500415

Overview of the benefits:

- ✓ Easy transfer between wheelchair - car - bed
- ✓ Ease of handling due to large grab handle
- ✓ Exceptionally hygienic due to a smooth gliding surface without grooves
- ✓ Easy to disinfect with all common disinfection agents
- ✓ Discreet and unobtrusive design
- ✓ Load bearing capacity up to 180 kg
- ✓ CE-certified medical device with aid number

Transa Sit On Towel

30,5 x 58 cm

Item no. 4410525

The Transa Sit On Towel supports the glide of the silver glider and simplifies the transfer with bare skin.

Anti-Glider Pads

M silver glider

14,5 x 18 cm

Item no. 4500418

L silver glider

20 x 20 cm

Item no. 4500419

The anti-glider pads are barely visible and very durable. They can be replaced as needed..

Beasy – the patient glide board with the turntable disc

In order to avoid shearing and friction forces, the Beasy models rest securely on the same spot, namely the transfer disc. The transfer disc always glides and rotates smoothly even under the maximum load.

Overview of the benefits:

- ✓ Easiest transfer due to gliding turntable disc
- ✓ Easy transfer between wheelchair - car - bed
- ✓ Ease of handling due to large grab handle on the rear side
- ✓ Clearly specified transfer path
- ✓ Particularly suitable for heavy, elderly and/or frail patients
- ✓ Considerable convenience for users, even with the help of another person
- ✓ Easy to disinfect with all common cleaning agents
- ✓ Load-bearing capacity up to 200 kg
- ✓ CE-certified medical device



Beasy Glider

80 x 26 cm
Item no. 4411080

Carrying Bag
Item no. 4411501

Beasy II

68 x 30 cm
Item no. 4411068

Carrying Bag
Item no. 4411500



Beasy Trans

102 x 30 cm
Item no. 4411102

Carrying Bag
Item no. 4411502

meinHOLTZ Transfer Board

65 x 24 cm



The meinHOLTZ Transfer Board - the design product with individual motifs

The meinHOLTZ is more than just a glide board. It combines the functionality and ease of use necessary for a transfer aid with appealing, modern designs. These traditionally rather plain aids are now available in great new designs. We were keen to develop a product that combines design and optimal functionality, and which is a pleasure to show to everyone who uses it.

Overview of the benefits:

- ✓ Modern designs from simple to colorful
- ✓ Easy transfer between wheelchair - car - bed
- ✓ A shape developed by therapists
- ✓ Minimized risk of slipping due to the gentle recess on the upper side
- ✓ Optimal stability due to multiple compressions
- ✓ Meets sustainability requirements due to wood as a renewable raw material
- ✓ Easy to disinfect with all common cleaning agents
- ✓ Load-bearing capacity up to 150 kg
- ✓ CE-certified medical device with aid number



Small DIY Tip:

Commercially available permanent markers can be exceptionally well applied to the glide board—For limitless creativity.

“The workload in the nursing professions is constantly growing. This also makes it all the more important to provide relief - for example, in the simple and ergonomic transfer and positioning of patients.”

Kai Stolpe, General Manager



Our promise of quality to you

User-safety is the top priority for us. As manufacturer and distributor of medical devices, we have been inspected and are certified according to DIN EN ISO 13485. All of our Transaroll medical devices are CE-marked according to the Medical Device Regulation (MDR). By placing our Rollboards, Transfer boards and further products on the market as medical devices, we guarantee that our products fulfill the essential safety requirements and the prescribed compliance procedures were performed. Thus, our products also fulfill the requirements of the Medical Device Law Implementation Act and the Medical Devices User Notification and Information Ordinance.



Transatlantic
Handels-Gesellschaft
Stolpe & Co. mbH

Siemensstraße 7
61267 Neu-Anspach

Tel: +49 (0) 60 81 94 30 50
Fax: +49 (0) 60 81 94 30 80

Email: info@transat.de
www.transatlantic.de