

# Transafix®

Fixation systems

## Medical Tubing – Simple, Hygienic & Securely Fixed

- ✓ Fixation strips
- ✓ Tube fixations
- ✓ Nasal tube fixations
- ✓ IV cannula fixation

 **Transatlantic**  
Products for a healthier world

# Made in Germany for more than 30 years

Our years of experience in the area of fixing systems is reflected in a particularly well-balanced and practice-oriented range of products. In addition to our fixation strips with a practical Velcro tab, we also offer nasal tube fixations for gastric or oxygen tubes, Endotracheal tube fixations for secure and hygienic intubation as well as classic catheter fixations.

## Fixation systems —the perfect solution for all concerns

Our fixations are found in the intensive care unit, in inpatient and home care, as well as the emergency services and are suitable for infants, children and adults. Patient comfort and simple, hygienic handling are our top priority. Therefore, our fixation systems are reviewed on a regular basis and continuously improved according to the latest technical advances. Since hygiene during patient care is of highest priority, we offer almost all of the products in our product range in sterile and non-sterile versions.



**Gentle even on sensitive skin**



**Maximum patient comfort**



**Suitable for use up to 30 days**



**Highest hygienic safety**

## About Transatlantic

Reliable cooperation is built on trust. For that reason, you need a partner you can trust. At Transatlantic, we have advised and supported our customers for over 85 years, and our mission is to create better individualized, more humane and efficient medical care for the benefit of patients and users. You will experience competent support and advice by our qualified medical product advisors when selecting products appropriate for you.



Transafix



Transafix



Transafix



Transafix



**Areas of application** 4




**Fixation strip** 5



**ET tube fixation** 8



**Nasal tube fixation** 9



**IV cannula fixation** 10



**Products - original size** 11

# Areas of application & products



## For intensive care units

The intensive care unit is the place where critically ill patients must be closely monitored and/or treated. This means that the nursing staff must devote their full attention to caring for the patient. The small details must happen simply, quickly and safely. For that reason, we have developed fixations which make everyday clinical routine easier for the nursing staff, providing a high level of patient comfort and safety. All product fixations for the intensive care unit are available sterile and thus meet the highest hygiene standard.

- ✓ Fixation strips
- ✓ ET tube fixation



## For inpatient care

There is often little time for patients during inpatient care in nursing homes or assisted living facilities as well as in hospital units. Therefore, it is important that processes are quick, safe and easy. Our fixations used during inpatient care are designed to be easy to handle and provide a high level of patient comfort.


- ✓ Fixation strips
- ✓ IV cannula fixation
- ✓ Nasal tube fixation



## For home care

In the home-care setting, patients often are not cared for by professional caregivers but by relatives. Therefore, they need products there that are easy to handle while being safe and hygienic.

- ✓ Fixation strips
- ✓ Nasal tube fixation



## For emergency services:

Emergency services carry out the initial care in emergency cases. Therefore, the procedures must flow in a rapid, hygienic and uncomplicated manner. Our products are ideally suited for this purpose since they are easy to handle and meet the highest hygiene and quality standards.

- ✓ Fixation strips
- ✓ ET tube fixation
- ✓ IV cannula fixation



## Transafix fixation strips – fixation for medical tubing with Velcro tab

Transfix fixation strips have been in use for more than 30 years in all areas of care. The basic carrier of the Velcro strap is the Transpore carrier patch, which is suitable for long-term applications\*. It includes an integrated Velcro fastener, with which the supply lines can be fixed to an adhesive surface.

### Transafix Fixation Strip No. 1

#### Sterile

50 pieces / box  
7.5 x 1.6 cm  
Item no. 8800100  
PPN: 111126903089

#### Non-sterile

100 pieces / box  
7.5 x 1.6 cm  
Item no. 8700150  
PPN: 110050611354

### Overview of the benefits:

- ✓ Freely movable or firm fixation of the tube in the Velcro tab
- ✓ Gentle on the skin, even sensitive skin
- ✓ Recommended, e.g., for PEG tubes, urinary catheters, CAPD, ...
- ✓ Suitable for long-term treatment\*
- ✓ Maximum patient comfort
- ✓ Highest quality standard
- ✓ Frequent opening and closing of the Velcro tab possible, as necessary



#### Areas of application

Smaller volumes, tubes, cables, lines



#### Service life

5 -7 days

\* long-term application > 30 days

## Transafix Fixation Strip No. 2

### Sterile

50 piece / box

9 x 3 cm

Item no. 8800200

PPN: 111126908264

### Non-sterile

50 piece / box

9 x 3 cm

Item no. 8700250

PPN: 110050610788



### Areas of application

PEG, CAPD, infusions, etc.



### Service life

5-7 days



### Areas of application

Bladder catheter/urinary catheter and suprapubic catheter as well as wound drainage, central venous catheter, lines, etc.



### Liegedauer

5-7 days

## Transafix Fixation Strip No. 3

### Sterile

30 pieces / box

15 x 4.5 cm

Item no. 8800300

PPN: 111126896442

### Non-sterile

25 pieces / box

15 x 4.5 cm

Item no. 8700350

PPN: 111126896442



## Transafix PEG Set No. 10

Item no. 8801004  
30 pieces / box

### Includes:

1 Transafix fixation strip No. 2  
1 drainage compress  
1 fixation pad 12 x 15 cm  
4 gauze compresses 7.5 x 7.5 cm  
PPN 11552871259

## Transafix PEG Set No. 10 – the complete solution for PEG tubes

With PEG tubes, there is a risk that they may be accidentally removed as a result of movements or of an excessive tensile load, which may result in inflammations at the puncture site. Therefore, they are often fixed with a simple suture and an additional tape patch. With our PEG Set No. 10 we have developed a complete solution that prevents this risk and is highly comfortable for the patient.

### Overview of the benefits:

- ✓ Complete solution, no additional products necessary
- ✓ Inflammations at the puncture site may be avoided by secure fixation of the PEG tube
- ✓ Long-term wearing period without skin irritations
- ✓ Safe and easy handling with maximum patient comfort




#### Areas of application

PEG tubes



#### Service life

5 –7 days (Transafix No. 2)



## Transafix ET Tube Fixation Set

Secure and comfortable fixation is important, especially during the application of Endotracheal tubes. Our Transafix ET Tube Fixation Set offers the perfect solution for this purpose. The set consists of two parts, 2 cheek pads, available in non-woven or hydrocolloid, as well as a soft, 40 cm long fixation strap. The holder can be individually adjusted and thus is suitable for all sizes.

The cheek pads and the fixation strap are connected by a secure Velcro fastener. Due to materials that are gentle on the skin, long-term application of up to 12 days is possible. Both materials, non-woven and hydrocolloid, are well-tolerated by the skin. We recommend using hydrocolloid cheek pads for patients with very sensitive or damaged skin.

### Overview of the benefits:

- ✓ Easy and secure fixation by the Velcro fastener on cheeks pads and fixation strap.
- ✓ Highest quality standard
- ✓ Gentle even on sensitive skin
- ✓ Quick and easy oral care possible
- ✓ More efficient than traditional solutions
- ✓ No residues on the tube
- ✓ Fixation possible without raising the head



#### Areas of application

Endotracheal tubes, feeding tubes, oxygen tubes, etc.



#### Service life

Non-woven c. 3 days, hydrocolloid c. 5 days, fixation strap as needed

### Transafix ET Tube Fixation – Non-woven set

Set with 2 cheek pads + 1 fixation strap

#### Sterile

Item no. 8800600  
30 pieces / box

#### Non-sterile

Item no. 8700600  
50 pieces / box  
PPN 110050609994

### Transafix ET Tube Fixation – Hydrocolloid set

Set with 2 cheek pads + 1 fixation strap

#### Sterile

Item no. 8600610  
30 pieces / box  
PPN 110948200634

### Transafix fixation tape

#### Sterile

Item no. 8800800  
50 pieces / box  
PPN 110948201200

#### Non-sterile

Item no. 8700800  
50 pieces / box  
PPN 11005060820

### Transafix non-woven cheek pads

#### Sterile

Item no. 8800700  
50 pairs / box

#### Non-sterile

Item no. 8700700  
50 pairs / box  
PPN 110050613610





## Transafix Nasal Tube Fixation

The use of nasogastric tubes is a frequent application in hospitals. The tube is inserted through the nose of the patient into the stomach. Therefore, the tube must be fixed securely on the bridge of the nose. Our Transafix Nasal Fixation System offers the perfect solution for this purpose. The combination of nasal holder, available in skin-tone non-woven or hydrocolloid, and the nasal strap make the system ready to use.

Both components are connected by a secure Velcro fastener. As a result, only the nasal strap needs to be replaced when the nasogastric tube is changed. Due to materials that are gentle on the skin, long-term application up to 12 days is possible. This reduces skin irritations and inflammations and prevents intensive follow-up care. Both non-woven and hydrocolloid are very well-tolerated by the skin. We recommend using a hydrocolloid nasal holder for very sensitive skin.

### Overview of the benefits:

- ✓ Easy connection and separation of the nasal holder and nasal strap by the Velcro fastener.
- ✓ Adaptable to the needs of the patient
- ✓ Easy removal without residues
- ✓ Gentle on the skin, even sensitive skin
- ✓ Highest quality standard
- ✓ Nasal holder is NOT removed during the tube change. Only the nasal strap is replaced.
- ✓ Tube is movable, preventing a pressure sore on the interior wall of the nose.



#### Areas of application

Feeding tubes, oxygen tubes, etc.



#### Service life

Non-woven c. 3 days, hydrocolloid approx. 5 days, nasal strap as needed

### Transafix Nasal Holder made of Hydrocolloid

#### Non-sterile

Item no. 8210100  
50 pieces / box

### Transafix Nasal Holder made of Non-woven

#### Non-sterile

Item no. 8210200  
50 pieces / box

### Transafix Nasal Strap

#### Non-sterile

Item no. 8210300  
100 pieces / box



#### IV cannula fixation

##### Sterile

100 pieces / box

Item no. 8900001P

PPN 111837942524

## IV Cannula Fixation

Our cannula fixation with non-adhesive wound dressing features outstanding adhesion properties and is very well-tolerated by the skin. As a result, it may remain on the skin for up to 7 days. It can be removed pain- and residue-free.

### Overview of the benefits:

- ✓ Easy, safe and comfortable fixation of indwelling peripheral cannulas
- ✓ Suitable for long-term use due to material being well-tolerated by the skin
- ✓ Easily removed without residues
- ✓ Ease of handling
- ✓ Maximum patient comfort
- ✓ Highest quality standard



#### Areas of application


Peripheral indwelling cannulas



#### Service life

5 -7 days

# Product illustrations - original size



*“Hygienic, safe, gentle on the skin.  
Simply more comfort for your patients  
without the extra cost.”*

**Kai Stolpe**, General Manager



## Our promise of quality to you

User-safety is the top priority for us. As manufacturer and distributor of medical devices, we have been inspected and are certified according to DIN EN ISO 9001 and DIN EN ISO 13485. All of our Transafix products are CE-marked according to the Medical Device Regulation (MDR) and have in part a Pharmacy Product Number (PPN). By placing our fixation systems on the market as medical devices, we guarantee that our products fulfill the essential safety requirements and the prescribed compliance procedures were performed. Thus, our products also fulfill the requirements of the Medical Device Law Implementation Act and the Medical Devices User Notification and Information Ordinance.



**Transatlantic  
Handelsgesellschaft  
Stolpe & Co. mbH**

Siemensstraße 21-23  
61267 Neu-Anspach

Tel: 06081 943050  
Fax: 06081 943080

Email: [info@transat.de](mailto:info@transat.de)  
[www.transatlantic.de](http://www.transatlantic.de)