

The syringes we offer are compatible with the following injectors – an overview

Whether CT or MRI – the Transaflow® Multi-Patient-Syringe-System for administration of contrast media enables a hygienic, safe, time-saving and resource-efficient working day in radiology. Both the syringes and the Multi-APS Systems are compatible with all standard syringes injectors.



Manufacturer	Device Type	Multi-APS & Patient Lines	Syringe	Syringe Specification		Application Time	Art.Number	PN Manufacturer
Medrad® Stellant	CT	✓	✓	Syringe	200 ml	up to 24 hrs up to 12 hrs	2850100 2840100	R100102 M100102
Medrad® Spectris	MRT	✓	✓	Syringe Set	2 x 65 ml	up to 24 hrs up to 12 hrs	2850001 2840001	M130103 M130102
Medrad® Spectris Solaris EP	MRT	✓	✓	Syringe Set	65/115 ml	up to 24 hrs up to 12 hrs	2850002 2840002	M130204 M100302
Medtron Accutron®	CT/MRT	✓	✓	Syringe	200 ml	up to 24 hrs up to 12 hrs	2850300 2840300	R400104 M400104
Nemoto Dual Shot alpha 7	CT	✓	✓	Syringe	100 ml	up to 24 hrs up to 12 hrs	2850411 2840411	R300101 M300101
Nemoto Dual Shot alpha 7	CT	✓	✓	Syringe	200 ml	up to 24 hrs up to 12 hrs	2850412 2840412	R300102 M300102
Nemoto SonicShot 7	MRT	✓	✓	Syringe	60 ml	up to 24 hrs up to 12 hrs	2850400 2840400	R340300 M340300
Guerbet Opti-Vantage® (before Mallinckrodt)	CT	✓	✓	Syringe	200 ml	up to 24 hrs up to 12 hrs	2850500 2840500	R200101 M200101
Bracco Empower CTA®	CT	✓	✓	Syringe Set	2 x 200 ml	Single-use Product	2820200	420102

Standardization as well as safety and hygiene in radiology routine


Our Multi-Patient-Syringe-System optimizes contrast medium administration routines. The system consists of our proven 12- or 24-hours Transaflow® Multi-APS-Contrast media-Filling System and the matching 12- or 24-hours syringes. The closed system is complemented by our PWL (straight) and PWLS (spiraled) patient lines in the desired length.

- ✓ With special drip and clog stop
- ✓ Multiple application up to 12/24 hrs
- ✓ Certified by the DGKH e.V.

For more information about our Transaflow® range, please contact: info@transat.de

Transatlantic – your expert for patient transfer, fixation systems, contrast media application and hospital supplies

Discover our product ranges for patient transfer, fixation systems and consumables for the daily hospital routine. Transatlantic offers products that are tailored to the specific requirements of clinics, care facilities, and medical supply stores.



Our promise to you: Products for a healthier world

Reliable cooperation is based on trust. At Transatlantic we have been advising and supporting our customers for over 85 years, and our mission is to create better, more individualized, more humane and more efficient medical care for the benefit of patients and users.

Benefit from our many years of experience in manufacturing and distributing our products. Our products are of the highest quality and are carefully selected and tested by us. We offer our customers a wide range of products tailored to the specific requirements and needs of their facilities. With our service, we want to help our customers make their work more effective and efficient, thereby contributing to the improvement of patient care.

Our qualified medical device consultants provide expert support and advice to help you select the right products.

- ✓ Personal consultation
- ✓ Customized solutions
- ✓ Excellent service
- ✓ RollboardCheck



Transatlantic
Handels-Gesellschaft
Stolpe & Co. mbH

Siemensstraße 7
61267 Neu-Anspach

Phone: 06081 943050
Fax: 06081 943080

Mail: info@transat.de
www.transatlantic.de