

Transaroll®

Patient Transfer

Practical Aids for Easy and Ergonomic Patient Transfer

- ✓ Roll boards & accessories
- ✓ Patient scale
- ✓ Storage systems
- ✓ Sheets
- ✓ Transfer boards & accessories

 **Transatlantic**
Products for a healthier world

Your expert for all types of transfer and positioning aids for more than 30 years!

In the 1980s, we introduced rollboards and sliding boards to the German market. Today, they are the standard of care for lateral and seated transfers in all clinical and care settings.

With easy handling, back-protecting working and optimal patient comfort, our Transaroll Boards and transfer boards facilitate everyday life in nursing and emergency services.

As a manufacturer located in Germany, we guarantee the highest quality and short delivery routes. In addition, we are DIN EN ISO 13485 certified. Our medical devices have undergone the relevant conformity assessment procedures and thus bear the CE mark as per the Medical Device Regulation (MDR).

Benefit from our service

To make it easier for you to implement the Medical Devices Operator Ordinance, we recommend our non-binding „Rollboard Check“ service. As part of a safety check, we suggest to inspect and identify the functionality and proper condition of your rollboards.

Training and instruction in proper handling are included. To make the work of the hospital staff even easier, they will receive a list of the rollboards with details of their condition and location. So they have the perfect overview of all of the various rollboards at their facility.



Protects the user's back



Cost-effective due to long product shelf life



Time-saving for more time for the patients



Over 30 years of experience with transfer and positioning aids



Extremely hygienic thanks to ease of cleaning or the use of special protective covers



Dignified procedure for the patient

About Transatlantic

Reliable cooperation is built on trust. For that reason, you need a partner you can trust. At Transatlantic, we have advised and supported our customers for over 85 years, and our mission is to create better individualized, more humane and efficient medical care for the benefit of patients and users. You will experience competent support and advice by our qualified medical product advisors when selecting products appropriate for you.



Rollboards & accessories **4**

Aids for lateral transfer



Patient scale **9**

Simply weigh patients in supine position



Storage systems **10**

Hygienic storage and easy transport of the rollboards



Sheets **11**

Absorbent transfer aids



Transfer boards & accessories **12**

Transfer aids for the wheelchair

Rollboards – aids for lateral transfer

	Silver	Hygiene	Blue	Professional	Smart	Klipp-Klapp
Sizes	5	3	4	5	2	1
Ends	Integrated carrying handles	Integrated carrying handles	Integrated carrying handles	Plastic end strip with loops	Plastic handle bar	none
Load	350 kg	200 kg	200 kg	200 kg	180 kg	150 kg
Radiolucent	✓	✓	✓	✓	✓	✓
Durability	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★	★ ★ ★	★ ★	★ ★
Gliding ability	★ ★ ★ ★ ★	★ ★ ★	★ ★ ★ ★	★ ★ ★ ★	★ ★ ★	★
Possible applications	★ ★ ★ ★ ★	★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★	★ ★ ★	★
Hygiene	★ ★	★ ★ ★ ★ ★	★ ★ ★	★ ★	★ ★	★ ★ ★ ★
	#1	#2	#3	#4	#5	#6

For more information on the available variants, please refer to the individual product information or our website at: www.transatlantic.de



What does this mean? Standard, Mini, Ambu, Long, Special and Slim are the size units of our rollboards.

**Please note that the sizes of all rollboards may vary by +/- 3%.*

Standard: ca. 180 x 50 cm.

Our most common size, since the boards can be used in all areas.

Mini: between 95-82 x 50 cm.

The boards are primarily used in gynecology and proctology, as they allow transfer even when the legs are elevated.

Ambu: between 86-75 x 45-40 cm.

Because of their narrow size, they are particularly suitable for short distances and transfers such as from narrow stretchers in emergency services. In addition, they are very space-saving when stored and can also be stored under a stretcher.

Long: approx. 195 x 50 cm.

Particularly suitable for large patients and for beds without disturbing the headboard and footboard.

Special: approx. 150 x 60 cm.

Particularly suitable for obese patients and transfers where a longer distance must be crossed.

Slim: approx. 120 x 40 cm.

Particularly suitable for stretcher transfers during emergency room admission and in CT.

For a smooth patient transfer, the heaviest points of the body—shoulders and pelvis—must be on the rollboard. That's why our Transaroll boards are most common in the Standard size. They can be used in all areas and allow for a comfortable transfer for both patient and user, since both the head and the heels rest on the board. But our portfolio also provides special solutions for the operating room with leg elevation (Mini), for the outpatient clinic and emergency room (Ambu, Slim), and for very tall or obese patients (Long, Special).

Basically, however, it is the width of the rollboard that determines its suitability for use. The wider the rollboard, the greater the distances that can be covered. Whether a rollboard is foldable or not has no influence on its transfer capabilities. They can only be stored in a more space-saving way, such as in our wall mount organizer.



Transaroll Silverboard - the rollboard with the silver cover

The Transaroll Silverboard is the new standard for all lateral positioning and transfer applications. Due to its extremely stable core, the Silverboard allows for ergonomic and easy patient transfer.

Standard (foldable)
180 x 50 cm
Item no. 4500030

Special (not foldable)
150 x 60 cm
Item no. 4500034

Mini (not foldable)
95 x 50 cm
Item no. 4500029

Special (foldable)
150 x 60 cm
Item no. 4500033

Mini (foldable)
95 x 50 cm
Item no. 4500028

Key facts

Core

Strong polyethylene core

Cover

Hygienic, antibacterial gliding fabric

Ends

Integrated carrying handles, extremely stable

Load


350 kg max.

Overview of the benefits:

- ✓ Ideal for obese patients
- ✓ The only rollboard which makes it possible to lift a person
- ✓ No sharp edges, all corners are rounded off
- ✓ No magnetic parts, plastic parts or screws
- ✓ Radiolucent and suitable for MRI
- ✓ Transfer even with soft positioning mattresses
- ✓ No end strip (maximum hygienic safety)
- ✓ Can be disinfected using all common disinfection agents

Transaroll Hygieneboard - maximum hygienic safety


The Hygieneboard ensures the most hygienic use of all. The rollboard consists of a stable and easy-to-clean core and our green gliding cover, which is applied to the rollboard each time it is used.




Standard (foldable)
180 x 50 cm
Item no. 4500026

Mini (foldable)
95 x 50 cm
Item no. 4500036

Mini (non-foldable)
95 x 50 cm
Item no. 4500035

Key facts	 Core	 Cover	 Ends	 Load
	Strong polyethylene core	Disposable gliding cover	Integrated carrying handles, extremely stable	200 kg max.



Transaroll Blueboard - the rollboard for care





Adapted to the average size and weight of patients in Europe, it is suitable for low-impact transfer in home care, as well as for use in clinics, hospitals and nursing homes that rely on quality and experience. The cover of the rollboard is machine washable.


Standard (foldable)
180 x 50 cm
Item no. 4501100

Long (foldable)
195 x 50 cm
Item no. 4500406

Mini (not foldable)
95 x 50 cm
Item no. 4500400

Mini (foldable)
95 x 50 cm
Item no. 4500403

Key facts	 Core	 Cover	 Ends	 Load
	Strong polyethylene core	Ripstop nylon	Integrated carrying handles, extremely stable	200 kg max.



Transaroll Professional board – the classic product for more than 30 years

Our Transaroll Professional boards have been on the German medical market for over 30 years and are as reliable today as they were then. They consist of a special lightweight foam, which is also very durable and resilient despite its low weight.

Standard (not foldable)
180 x 50 cm
Item no. 4401150

Standard (foldable)
180 x 50 cm
Item no. 4401100

Mini (not foldable)
82 x 50 cm
Item no. 4400400

Ambu (not foldable)
80 x 40 cm
Item no. 4400114

Slim (not foldable)
120 x 40 cm
Item no. 44004200

Key facts



Core

Special lightweight foam



Cover

Hygienic, antibacterial gliding sheet



Ends

Plastic endstrip with loops



Load

200 kg max.

Transaroll Smartboard - the low-cost alternative

This rollboard is perfect for those who are looking for an inexpensive rollboard, but still want to enjoy high gliding capability.



Standard (foldable)
180 x 50 cm
Item no. 4400102

Ambu (foldable)
86 x 40 cm
Item no. 4400101

Key facts



Core

Plate welded in polyester fabric



Cover

Hybrid nylon cover



Ends


Plastic handle bar



Load





180 kg max.

Transaroll Klipp-Klapp Board - the roll board for every ambulance



The Transaroll Klipp-Klapp Board was developed with emergency services in mind. Due to the double foldable body of the board, it allows for quick and easy space-saving storage in the ambulance.

Ambu (double foldable)
75 x 45 cm
Item no. 4400103

Key facts	 Core HDPE Board	 Cover Grey gliding fabric	 Ends No grip bars	 Load max. 150 kg
------------------	---	---	--	--

Accessories - Covers for fast, hygienic use



Disposable Gliding Cover

Due to the exceptional gliding properties, the cover is suitable for moving the patient in bed even without a rollboard.

Standard, 60 pcs.
(Roll, perforated)
180 x 53 cm
Item no. 4400204

Disposable, by the meter
Roll, can be cut to size)
100 m x 60 cm
Item no. 4400200

Disposable, by the meter
(cardboard)
100 m x 60 cm
Item no. 4400203




Disposable Protective Cover

Protects the rollboard from possible contamination in individual cases.

Disposable protective cover, 150 pc, anti-static
(Roll, perforated)
54 x 100 cm
Item no. 4500126

Disposable protective cover, 100 pc, anti-static
(Roll, perforated)
56 cm x 195 cm
Item no. 4500125N

Disposable protective cover, 50 pc, standard
(Cardboard, folded)
60 cm x 195 cm
Item no. 4500120



Patient scale - simply weigh lying down

With the innovative mobile patient scale, patients can be easily weighed while lying down or during transfer. Thanks to the fast and efficient weighing process, which can be easily integrated into daily workflows, weight no longer needs to be estimated. Accurate weight information helps to dose medications precisely and thus also save lives.

180 x 70 x 3 cm
Item no.: 4500900

Key facts

- | | |
|--|--|
| <p> Calibration class
Calibration class 3</p> | <p> Cover
PVC-coated polyester mat</p> |
| <p> Grips
Integrated carrying handles</p> | <p> Load
250 kg max.</p> |

Overview of the benefits:

- ✓ Transfer board with built-in scale
- ✓ Weighs the patient quickly and accurately during transfer
- ✓ Dignified weighing, less stress for the patient
- ✓ Quick, easy, time-saving
- ✓ Accurate administration of weight-dependent medication
- ✓ Large display with hold function
- ✓ Capacity: 250 kg
- ✓ Rechargeable battery with 50 hours runtime during continuous operation

Storage systems - hygienic storage and ease of handling of rollboards and more

With our practical storage systems, all the necessary products for a quick, reliable and hygienic transfer and positioning are always at hand. The rollboards do not touch the floor when the disposable covers are applied, and cross-contamination can be avoided.



Organizer – storage for rollboard and cover

The organizer can be easily attached to any door or wall. The upper compartment of the organizer has an outlet for the disposable covers, while the lower compartment can be used to conveniently store the rollboards on edge.

B 70,5 x H 119 x T 17 cm
Item no. 4500417

Transfer Trolley – storage for patient scale, rollboard & cover

The Transfer Trolley not only provides storage options for the patient scales, the rollboard and the disposable covers but also serves as a practical transport device to quickly bring the transfer and positioning aids to their place of use. Without any heavy lifting.

Item no. 4500416



Sheets - absorbent transfer aids


The revolutionary sheets are not only used for lifting, carrying and extricating patients, but are also ideal as waterproof padding and for transferring a patient in a way that protects the back.

Overview of the benefits:

- ✓ Padding and transfer cloth with handles in one product
- ✓ Safe and stable gripping due to multiple handles
- ✓ Can be used alone or together with our Transaroll boards
- ✓ Quick and clean use
- ✓ Disposable product


Transatlantic Multisheet

80 x 210 cm
Item no. 4700181



Key facts


<p> Material</p> <p>Outer surface: Polyethylene cover, blue, breathable</p> <p>Absorbent core: Superabsorbent, fluff- tissue mixture</p> <p>Fleece: Polypropylene (water absorbent or water permeable)</p>	<p> Grip holes</p> <p>8 integrated grip holes</p> <p> Fluid uptake</p> <p>2.8 l fluid uptake</p> <p> Load</p> <p>150 kg max.</p>
--	--



XXL Carry Sheet

107 x 210 cm
Item no. 4700081

Key facts

<p> Material</p> <p>Polyester antimicrobial fabric</p>	<p> Fluid uptake</p> <p>4 l fluid uptake</p>
<p> Grip holes</p> <p>16 integrated grips</p>	<p> Load</p> <p>450 kg max.</p>

Transfer boards & accessories - transfer aids for the wheelchair

Our broad portfolio also offers useful aids for seated transfers. Our glide boards are ideally suited for all care areas in clinical and nursing care applications and facilitate ergonomic working.




Silver glider

Classic glider for all applications



Beasy glide board

Nursing glide board with turntable disc



meinHOLTZ transfer board

Design product with individual motifs

What is a transfer board?

Transfer boards make it easier to transfer mobility-impaired people. They form a bridge between the wheelchair and another surface, such as the bed, toilet or car seat.

Who can use a transfer board?

In principle, anyone who finds it difficult to transfer from a wheelchair can use a glide board. It supports an active transfer with the greatest possible independence. A glide board can also be a great support for relatives or people who help people in need of care, and can relieve strain on the back.



Silver Glider - the classic glider in all areas

The M and L Silver Gliders have a shape developed by therapists, which has proven its worth for many years. Due to a very smooth surface, the gliders minimize frictional resistance while anti-slip pads prevent the board from slipping away on the base.



M Silver Glider

67 x 22 cm

Item no. 4500414

L Silver Glider

85 x 25 cm

Item no. 4500415

Overview of the benefits:

- ✓ Easy transfer between wheelchair - car - bed
- ✓ Ease of handling due to large grab handle
- ✓ Exceptionally hygienic due to a smooth gliding surface without grooves
- ✓ Easy to disinfect with all common disinfection agents
- ✓ Discreet and unobtrusive design
- ✓ Load bearing capacity up to 180 kg
- ✓ CE-certified medical device with aid number



Transa Sit On Towel

30,5 x 58 cm

Item no. 4410525

The Transa Sit On Towel supports the glide of the silver glider and simplifies the transfer with bare skin.

Anti-Glider Pads

M silver glider

14,5 x 18 cm

Item no. 4500418

L silver glider

20 x 20 cm

Item no. 4500419


The anti-glider pads are barely visible and very durable. They can be replaced as needed..

Beasy – the patient glide board with the turntable disc

In order to avoid shearing and friction forces, the Beasy models rest securely on the same spot, namely the transfer disc. The transfer disc always glides and rotates smoothly even under the maximum load.

Overview of the benefits:

- ✓ Easiest transfer due to gliding turntable disc
- ✓ Easy transfer between wheelchair - car - bed
- ✓ Ease of handling due to large grab handle on the rear side
- ✓ Clearly specified transfer path
- ✓ Particularly suitable for heavy, elderly and/or frail patients
- ✓ Considerable convenience for users, even with the help of another person
- ✓ Easy to disinfect with all common cleaning agents
- ✓ Load-bearing capacity up to 200 kg
- ✓ CE-certified medical device



Beasy Glider

80 x 26 cm
Item no. 4411080

Carrying Bag
Item no. 4411501

Beasy II

68 x 30 cm
Item no. 4411068


Carrying Bag
Item no. 4411500



Beasy Trans

102 x 30 cm
Item no. 4411102

Carrying Bag
Item no. 4411502



meinHOLTZ Transfer Board

65 x 24 cm



The meinHOLTZ Transfer Board - the design product with individual motifs

The meinHOLTZ is more than just a glide board. It combines the functionality and ease of use necessary for a transfer aid with appealing, modern designs. These traditionally rather plain aids are now available in great new designs. We were keen to develop a product that combines design and optimal functionality, and which is a pleasure to show to everyone who uses it.

Overview of the benefits:

- ✓ Modern designs from simple to colorful
- ✓ Easy transfer between wheelchair - car - bed
- ✓ A shape developed by therapists
- ✓ Minimized risk of slipping due to the gentle recess on the upper side
- ✓ Optimal stability due to multiple compressions
- ✓ Meets sustainability requirements due to wood as a renewable raw material
- ✓ Easy to disinfect with all common cleaning agents
- ✓ Load-bearing capacity up to 150 kg
- ✓ CE-certified medical device with aid number



Small DIY Tip:

Commercially available permanent markers can be exceptionally well applied to the glide board—For limitless creativity.

“The workload in the nursing professions is constantly growing. This also makes it all the more important to provide relief - for example, in the simple and ergonomic transfer and positioning of patients.”

Kai Stolpe, General Manager



Our promise of quality to you

User-safety is the top priority for us. As manufacturer and distributor of medical devices, we have been inspected and are certified according to DIN EN ISO 9001 and DIN EN ISO 13485. All of our Transafix products are CE-marked according to the Medical Device Regulation (MDR) and have in part a pharmaceutical registration number (PZN). By placing our fixation systems on the market as medical devices, we guarantee that our products fulfill the essential safety requirements and the prescribed compliance procedures were performed. Thus, our products also fulfill the requirements of the Medical Device Law Implementation Act and the Medical Devices User Notification and Information Ordinance.



Transatlantic
Handelsgesellschaft
Stolpe & Co. mbH

Siemensstraße 21-23
61267 Neu-Anspach

Tel: 06081 943050
Fax: 06081 943080

Email: info@transat.de
www.transatlantic.de