

Transafix®

Fixation systems

Medical Tubing – Simple, Hygienic & Securely Fixed

- ✓ Fixation strips
- ✓ Tube fixations
- ✓ Nasal tube fixations
- ✓ IV cannula fixation



Made in Germany for more than 30 years

Our years of experience in the area of fixing systems is reflected in a particularly well-balanced and practice-oriented range of products. In addition to our fixation strips with a practical Velcro tab, we also offer nasal tube fixations for gastric or oxygen tubes, Endotracheal tube fixations for secure and hygienic intubation as well as classic catheter fixations.

Fixation systems —the perfect solution for all concerns

Our fixations are found in the intensive care unit, in inpatient and home care, as well as the emergency services and are suitable for infants, children and adults. Patient comfort and simple, hygienic handling are our top priority. Therefore, our fixation systems are reviewed on a regular basis and continuously improved according to the latest technical advances. Since hygiene during patient care is of highest priority, we offer almost all of the products in our product range in sterile and non-sterile versions.



Gentle even on sensitive skin



Maximum patient comfort



Suitable for use up to 30 days



Highest hygienic safety

About Transatlantic

Reliable cooperation is built on trust. For that reason, you need a partner you can trust. At Transatlantic, we have advised and supported our customers for over 85 years, and our mission is to create better individualized, more humane and efficient medical care for the benefit of patients and users. You will experience competent support and advice by our qualified medical product advisors when selecting products appropriate for you.



Areas of application **4**



Fixation strip **5**




ET tube fixation **8**



Nasal tube fixation **9**



IV cannula fixation **10**



Products - original size **11**

Areas of application & products



For intensive care units

The intensive care unit is the place where critically ill patients must be closely monitored and/or treated. This means that the nursing staff must devote their full attention to caring for the patient. The small details must happen simply, quickly and safely. For that reason, we have developed fixations which make everyday clinical routine easier for the nursing staff, providing a high level of patient comfort and safety. All products for the intensive care unit meet the highest hygiene standards, and some are also available in sterile versions.

- ✓ Fixation strips
- ✓ ET tube fixation
- ✓ Multi-Fixation



For inpatient care

There is often little time for patients during inpatient care in nursing homes or assisted living facilities as well as in hospital units. Therefore, it is important that processes are quick, safe and easy. Our fixations used during inpatient care are designed to be easy to handle and provide a high level of patient comfort.

- ✓ Fixation strips
- ✓ IV cannula fixation
- ✓ Nasal tube fixation
- ✓ Multi-Fixation



For home care

In the home-care setting, patients often are not cared for by professional caregivers but by relatives. Therefore, they need products there that are easy to handle while being safe and hygienic.

- ✓ Fixation strips
- ✓ Nasal tube fixation
- ✓ Multi-Fixation



For emergency services:

Emergency services carry out the initial care in emergency cases. Therefore, the procedures must flow in a rapid, hygienic and uncomplicated manner. Our products are ideally suited for this purpose since they are easy to handle and meet the highest hygiene and quality standards.

- ✓ Fixation strips
- ✓ ET tube fixation
- ✓ IV cannula fixation



Transafix fixation strips – fixation for medical tubing with Velcro tab

Transfix fixation strips have been in use for more than 30 years in all areas of care. The basic carrier of the Velcro strap is the Transpore carrier patch, which is suitable for long-term applications*. It includes an integrated Velcro fastener, with which the supply lines can be fixed to an adhesive surface.

Transafix Fixation Strip No. 1

Sterile

50 pieces / box
7.5 x 1.6 cm
Item no. 8800100
PPN: 111126903089

Non-sterile

100 pieces / box
7.5 x 1.6 cm
Item no. 8700150
PPN: 110050611354

Overview of the benefits:

- ✓ Freely movable or firm fixation of the tube in the Velcro tab
- ✓ Gentle on the skin, even sensitive skin
- ✓ Recommended, e.g., for PEG tubes, urinary catheters, CAPD, ...
- ✓ Suitable for long-term treatment*
- ✓ Maximum patient comfort
- ✓ Highest quality standard
- ✓ Frequent opening and closing of the Velcro tab possible, as necessary



Areas of application

Smaller volumes, tubes,
cables, lines



Service life

5 -7 days

* long-term application > 30 days

Transafix Fixation Strip No. 2

Sterile

50 piece / box

9 x 3 cm

Item no. 8800200

PPN: 111126908264

Non-sterile

50 piece / box

9 x 3 cm

Item no. 8700250

PPN: 110050610788




Areas of application

PEG, CAPD, infusions, etc.



Service life

5 -7 days



Areas of application

Bladder catheter/urinary catheter and suprapubic catheter as well as wound drainage, central venous catheter, lines, etc.



Service life

5 -7 days

Transafix Fixation Strip No. 3

Sterile

30 pieces / box

15 x 4.5 cm

Item no. 8800300

PPN: 111126896442

Non-sterile

25 pieces / box

15 x 4.5 cm

Item no. 8700350

PPN: 111126896442



Transafix PEG Set No. 10

Item no. 8801004
30 pieces / box

Includes:

1 Transafix fixation strip No. 2
1 drainage compress 6 x 8 cm
1 fixation pad 12 x 15 cm
4 gauze compresses 7.5 x 7.5 cm
PPN 11552871259

Transafix PEG Set No. 10 – the complete solution for PEG tubes

With PEG tubes, there is a risk that they may be accidentally removed as a result of movements or of an excessive tensile load, which may result in inflammations at the puncture site. Therefore, they are often fixed with a simple suture and an additional tape patch. With our PEG Set No. 10 we have developed a complete solution that prevents this risk and is highly comfortable for the patient.

Overview of the benefits:

- ✓ Complete solution, no additional products necessary
- ✓ Inflammations at the puncture site may be avoided by secure fixation of the PEG tube
- ✓ Long-term wearing period without skin irritations
- ✓ Safe and easy handling with maximum patient comfort



Areas of application

PEG tubes



Service life

5 –7 days (Transafix No. 2)



Transafix ET Tube Fixation Set

Secure and comfortable fixation is important, especially during the application of Endotracheal tubes. Our Transafix ET Tube Fixation Set offers the perfect solution for this purpose. The set consists of two parts, 2 cheek pads, available in non-woven or hydrocolloid, as well as a soft, 40 cm long fixation strap. The holder can be individually adjusted and thus is suitable for all sizes.

The cheek pads and the fixation strap are connected by a secure Velcro fastener. Due to materials that are gentle on the skin, long-term application of up to 12 days is possible. Both materials, non-woven and hydrocolloid, are well-tolerated by the skin. We recommend using hydrocolloid cheek pads for patients with very sensitive or damaged skin.

Overview of the benefits:

- ✓ Easy and secure fixation by the Velcro fastener on cheeks pads and fixation strap.
- ✓ Highest quality standard
- ✓ Gentle even on sensitive skin
- ✓ Quick and easy oral care possible
- ✓ More efficient than traditional solutions
- ✓ No residues on the tube
- ✓ Fixation possible without raising the head



Areas of application

Endotracheal tubes, feeding tubes, oxygen tubes, etc.



Service life

Non-woven c. 3 days, hydrocolloid c. 5 days, fixation strap as needed

Transafix ET Tube Fixation – Non-woven set Set with 2 cheek pads + 1 fixation strap

Sterile

Item no. 8800600
30 pieces / box
PNN 111881714475

Non-sterile

Item no. 8700600
50 pieces / box
PPN 110050609994

Transafix fixation tape

Sterile

Item no. 8800800
50 pieces / box
PPN 110948201200

Non-sterile

Item no. 8700800
50 pieces / box
PPN 110050608207

Non-woven cheek pads

Sterile

Item no. 8800700
50 pairs / box
PNN 111881715041

Non-sterile

Item no. 8700700
50 pairs / box
PPN 110050613610

Cheek pads – Hydrocolloid

Non-sterile

Item no. 8800701
50 pairs / box



Transafix Nasal Tube Fixation

The use of nasogastric tubes is a frequent application in hospitals. The tube is inserted through the nose of the patient into the stomach. Therefore, the tube must be fixed securely on the bridge of the nose. Our Transafix Nasal Fixation System offers the perfect solution for this purpose. The combination of hydrocolloid nasal holder and the nasal strap make the system ready to use.

Both components are connected by a secure Velcro fastener. Due to materials that are gentle on the skin, long-term application up to 12 days is possible. This reduces skin irritations and inflammations and prevents intensive follow-up care. Hydrocolloid is very well tolerated by the skin and is also recommended for sensitive skin.

Overview of the benefits:

- ✓ Easy connection and separation of the nasal holder and nasal strap by the Velcro fastener.
- ✓ Adaptable to the needs of the patient
- ✓ Easy removal without residues
- ✓ Gentle on the skin, even sensitive skin
- ✓ Highest quality standard
- ✓ Nasal holder is NOT removed during the tube change. Only the nasal strap is replaced.
- ✓ Tube is movable, preventing a pressure sore on the interior wall of the nose.



Areas of application

Feeding tubes, oxygen tubes, etc.



Service life

hydrocolloid approx. 5 days,
nasal strap as needed


Transafix Nasal Holder made of Hydrocolloid

Non-sterile

Set with 1 Nasal Holder + 1 Nasal -strap

Item no. 8210400

50 pieces / box



IV Cannula Fixation

The cannula fixation from the manufacturer Pharmaplast S.A.E. with non-adhesive wound dressing features outstanding adhesion properties and is very well-tolerated by the skin.

As a result, it may remain on the skin for up to 7 days. It can be removed pain- and residue-free.

IV cannula fixation Sterile

100 pieces / box

Item no. 8900001P

PPN 111837942524

Overview of the benefits:

- ✓ Easy, safe and comfortable fixation of indwelling peripheral cannulas
- ✓ Suitable for long-term use due to material being well-tolerated by the skin
- ✓ Easily removed without residues
- ✓ Ease of handling
- ✓ Maximum patient comfort
- ✓ Highest quality standard



Areas of application

Peripheral indwelling cannulas



Service life

5-7 days

Transafix Multi-Fixation

The Transafix Multifixation is an innovative solution for simple and hygienic fixation of medical tubing on beds, bed linen or clothing. The product consists of a robust adhesive plaster that has been specially developed to adhere to surfaces such as fabric or furniture rather than to the skin. There is a strong Velcro strip on the upper side that enables the secure attachment of all kinds of hoses. 20 cm long and 2.5 cm wide, the Transafix Multifixation offers sufficient space for organized hose and cable management.

Multi-Fixation Non-sterile

Item no. 8100200






25 pieces / box

Overview of the benefits:

- ✓ Tubing can be attached or detached quickly
- ✓ Minimizes contamination, keeps tubing clean and tidy
- ✓ For various types of tubing, ideal in hospitals, care homes or at home
- ✓ Easy to remove without leaving residue or damage
- ✓ Simple design for intuitive handling
- ✓ No metallic parts, therefore suitable for MRI



Product illustrations - original size




*“Hygienic, safe, gentle on the skin.
Simply more comfort for your patients
without the extra cost.”*

Kai Stolpe, General Manager



Our promise of quality to you

User-safety is the top priority for us. As manufacturer and distributor of medical devices, we have been inspected and are certified according to DIN EN ISO 13485. All of our Transafix products are CE-marked according to the Medical Device Regulation (MDR) and have in part a Pharmacy Product Number (PPN). By placing our fixation systems on the market as medical devices, we guarantee that our products fulfill the essential safety requirements and the prescribed compliance procedures were performed. Thus, our products also fulfill the requirements of the Medical Device Law Implementation Act and the Medical Devices User Notification and Information Ordinance.



**Transatlantic
Handels-Gesellschaft
Stolpe & Co. mbH**

Siemensstraße 7
61267 Neu-Anspach

Tel: +49 (0) 60 81 94 30 50
Fax: +49 (0) 60 81 94 30 80

Email: info@transat.de
www.transatlantic.de